


H-STERISCOPE™ I

SINGLE-USE BRONCHOSCOPES



ZERO

SLIM

NORMAL

LARGE

EXTRA



LARGE RANGE OF STERILE SINGLE-USE BRONCHOSCOPES

- No cross-contamination
- No repair & reprocessing cost



EASY INSERTION AND PROGRESSION

- Up and down bending up to 210°
- Left of right rotation up to 90°



OPTIMAL VISUALIZATION

- LED lighting
- 110° field of view

References	Model	Color	Working channel diameter (mm)	Insertion tube outer diameter (mm)
BCV1-O2*	Zero		No working channel	2.5
BCV1-C2	Slim		1.2	3.2
BCV1-M2	Normal		2.2	4.9
BCV1-S2	Large		2.8	5.8
BCV1-W2	Extra		3.2	6.2

* Available on special request.

ALL-IN-ONE VIDEO MONITOR



REFERENCE: DVM-A1



ALL IN ONE HIGH RESOLUTION TOUCH SCREEN

- Screen size: 12.1"
- Maximum monitor resolution: 1280x800 pixels
- Ergonomic and instinctive navigation



DATA RECORDS

- Storage capacity: 64 GB
- Saving via 2 USB ports
- Video format: MP4 / Photo format: png



POWER SUPPLY & CONNECTIONS

- Main power cord
- Lithium battery with 6 hours autonomy
- HDMI and DVI cords

